

# WILAméd

*Equipment for Professionals*



## VENTURA TFM Total Face Mask

Disposable patient interfaces  
for noninvasive ventilation

[www.wilamed.com](http://www.wilamed.com)



# VENTURA Total Full Face

WILMed's Total Face Masks have been designed and created with the patient's comfort in mind. Our VENTURA interfaces are less intimidating than other patient interfaces, as our Total Face Masks promote an unrivaled openness, enabling better treatment outcome with anxious claustrophobic patients.

## Product features at a glance:

- Single Patient Use (For up to 14 days)
- Three different (Elbow) versions, including headgear
- Available in three sizes (S, M, L)
- Accommodating mouth breathing and lack of teeth
- For use with all major noninvasive ventilators

## The **key benefits** for patient and care givers are numerous:

- Provides clear and unobstructed visibility in all directions
- Ultra soft surgical silicone cushion for a low pressure, soft fit and seal all around the face
- Advanced design to minimize facial pressure points during long duration of mechanical ventilation
- 4 point mask strap that snugly fits against the head's natural curvature
- No more annoying pressure on the forehead or nose

## The **additional benefits** are:

- Usable for up to 14 days (single patient use)
- Free from DEHP, PVC and Latex
- MRI Safe for use with MRI compatible ventilators
- For adult patients >20kg
- Wide range of critical and home care breathing circuits



**NV/ Standard Blue Elbow**



**Nonvented with AAV**



**Vented with AAV**

# VENTURA Total Full Face -

## Key benefits for patients and care givers

- Pure and clear mask body
- Soft Silicon mask cushion
- Anti asphyxia valve
- Quick Release mask tube connector
- ETCO2 Sampling port
- Nasogastric port (NG)
- Quick Release strap buckles
- 360° rotating mask elbow

Size	Weight (g)	Dead Space (ml)
L	189	535
M	180	483
S	149	315

## Order Information

### SBE Standard Blue Elbow (nonvented)

P/N	Description	Qty.
201112	VENTURA SBE Total Face, incl. Strap, Size S	1
201111	VENTURA SBE Total Face, incl. Strap, Size M	1
201110	VENTURA SBE Total Face, incl. Strap, Size L	1



### VENTED/AAV

P/N	Description	Qty.
201113	VENTURA Vented /AAV Total Face, incl. Strap, Size S	1
201114	VENTURA Vented /AAV Total Face, incl. Strap, Size M	1
201115	VENTURA Vented /AAV Total Face, incl. Strap, Size L	1



### NV/AAV ( nonvented with anti asphyxia)

P/N	Description	Qty.
201107	VENTURA NV/AAV Total Face, incl. Strap, Size S	1
201108	VENTURA NV/AAV Total Face, incl. Strap, Size M	1
201109	VENTURA NV/AAV Total Face, incl. Strap, Size L	1





## VENTURA TFM

### Total Face Mask

## WILAmed GmbH

Gewerbepark Barthelmesaurach  
Aurachhöhe 5–7  
91126 Kammerstein (Germany)

Phone: +49 9178 996999-0 Fax:  
+49 9178 996778  
[info@wilamed.com](mailto:info@wilamed.com)  
[www.wilamed.com](http://www.wilamed.com)

