


WILAméd

Equipment for Professionals



**WILAhiflow
Maxi**

www.wilamed.com

NHFOT Nasal Cannula for Adults



WILAhiflow Maxi

The **WILAhiflow Maxi** Nasal Cannula is intended to provide nasal high flow oxygen therapy (NHFOT) for adult patients with acute respiratory failure. The Nasal Cannula is designed to deliver air, with or without blended oxygen, up to a maximum flow of 60 l/min. Special feature of the cannula – the flexible smooth bore connecting tube. WILAhiflow Maxi Nasal Cannula is for single-patient use in the hospital/clinical and homecare environment.

Order Details

| Part Number | Item | Size | Box/Qty. |
|-------------|-------------------------|------|----------|
| 30076310 | WILAhiflow Maxi, small | S | 10 |
| 30076410 | WILAhiflow Maxi, medium | M | 10 |
| 30076510 | WILAhiflow Maxi, large | L | 10 |



Distributed by:

WILamed GmbH

Medizinische Geräte und Zubehör

Gewerbepark Barthelmesaurach
Aurachhöhe 5–7
91126 Kammerstein (Germany)

Phone: +49 9178 996999-0
Fax: +49 9178 996778
info@wilamed.com
www.wilamed.com

