


# WILAméd

*Equipment for Professionals*



## WILAhiflow Mini

[www.wilamed.com](http://www.wilamed.com)

NHFOT Nasal Cannula for Infants



# WILAhiflow Mini

The **WILAhiflow Mini** Nasal Cannula for Infant is a single use nasal cannula intended for use with a nasal high flow oxygen therapy (NHFOT) system to deliver heated and humidified nasal high flow therapy to spontaneously breathing patients who require breathing support. WILAhiflow Mini Nasal Cannula is designed for use in hospital environments and must be prescribed by a physician.

## Order Details

Part Number	Item	Color	Size	Box/Qty.
30076010	WILAhiflow Mini, small	Yellow	S	10
30076110	WILAhiflow Mini, medium	Purple	M	10
30076210	WILAhiflow Mini, large	Green	L	10



### Distributed by:

## WILAmEd GmbH

Medizinische Geräte und Zubehör

Gewerbepark Barthelmesaurach  
Aurachhöhe 5–7  
91126 Kammerstein (Germany)

Phone: +49 9178 996999-0  
Fax: +49 9178 996778  
info@wilamed.com  
www.wilamed.com

