


WILAméd

Equipment for Professionals



**WILAhiflow Mini
Plus Nasal Cannula**
www.wilamed.com



WILAhiflow Mini Plus Nasal Cannula


The **WILAhiflow Mini Plus** Nasal Cannula for Infant is a single use nasal cannula intended for use with a nasal high flow oxygen therapy (NHFOT) system to deliver heated and humidified nasal high flow therapy to spontaneously breathing patients who require breathing support. WILAhiflow Mini Nasal Cannula is designed for use in hospital environments and must be prescribed by a physician.

Adjustable Slider

Thanks to the slider, secure and comfortable fitting can be achieved without need of an additional strap.

Skin-friendly Stickers

The bilateral stickers are designed under the conception of stability, good air permeability and skin-friendly.



Soft-touching Tubing

Gentle touching on the ears and face so that the friction and pressure on delicate tissues are minimized.

Multiple Patient Groups

Color-coded cannula sizes fulfill requirements from various patient groups.

Order Details

Part Number	Item	Color	Flow Range	Patient's Weight Range
3007671	WILAhiflow Mini Nasal Cannula, small	Yellow	2 ~ 10 L/min	2 ~ 8 Kg
3007681	WILAhiflow Mini Nasal Cannula, medium	Blue	2 ~ 20 L/min	3 ~ 15 Kg
3007691	WILAhiflow Mini Nasal Cannula, large	Green	2 ~ 25 L/min	12 ~ 22 Kg

Distributed by:

WILAmEd GmbH

A company of the Löwenstein Group

Gewerbepark Barthelmesaurach
Aurachhöhe 5–7
91126 Kammerstein (Germany)

Phone: +49 9178 996999-0
Fax: +49 9178 996778
info@wilamed.com
www.wilamed.com



BMC Medical Co., Ltd.
Room 110, Tower A Fengyu Building,
No.115 Fucheng Road, Haidian
100036 Beijing, China

Tel: +86-22-8293 9881
Fax: +86-22-8293 9881 etx. 810
en.bmc-medical.com