

WILAméd

Equipment for Professionals



Product Catalogue 2024

As a Health Care Professional, you need information that helps you to select the right medical equipment for a wide variety of respiratory applications. This product catalogue will help you to choose from industry leading equipment for critical- and home care applications, making it easier to find the right product for the job.

Our extensive industry knowledge drives us toward developing answers to the medical industry challenges of today and tomorrow.



COMMITMENT TO QUALITY

Our products are built in world-class production facilities and are subject to rigorous testing exceeding industry standards for quality and performance.



SUSTAINABLE PRODUCTION

All WILAméd products meet relevant ISO and EN Standards, such as ISO 13485 and Annex II of the Medical Device Directive, which allows CE marking of all our products.

■ made
■ in
■ Germany




ALWAYS THE RIGHT SOLUTION

WILAméd - advanced respiratory care products, offering innovation, quality and safety.

WILAméd

Wir erleichtern das Atmen.



WILAméd is an internationally operating manufacturer and distributor of medical devices and accessories in the field of invasive and non-invasive ventilation. In over 60 countries around the world, WILAméd products make an important contribution to respiratory patient care for neonates, children and adults.

The medium-sized company serves the homecare sector as well as clinics, hospitals and rehabilitation facilities through authorized specialist trade partners.

In the Barthelmesaurach industrial park, more than 150 employees in production and administration work every day to optimize the treatment success of ventilated patients.

In the fight against the Corona pandemic, WILAméd is one of the system-relevant companies in the healthcare sector as part of Germany's critical infrastructure.

WILAméd GmbH

A Company of the
Löwenstein Group

Aurachhöhe 5–7
91126 Kammerstein
info@wilamed.com
www.wilamed.com

TABLE OF CONTENTS

ACTIVE RESPIRATORY GAS HUMIDIFIERS 2

DISPOSABLE DOUBLE LIMB BREATHING SETS - HEATED 6

DISPOSABLE NEONATAL HFO BREATHING SETS - HEATED..... 8

DISPOSABLE SINGLE LIMB BREATHING SETS - HEATED11

DISPOSABLE BREATHING SETS - NON-HEATED 14

HUMIDIFICATION CHAMBERS..... 15

DISPOSABLE BREATHING SETS HOME CARE - HEATED 17

DISPOSABLE BREATHING SETS HOME CARE - NON-HEATED 24

DISPOSABLE BREATHING TUBES HOME CARE - NON-HEATED..... 29

DISPOSABLE PATIENT INTERFACES..... 30

HEAD GEAR FOR PATIENT INTERFACES - HOME-CARE 33

MECHANICAL AIR / OXYGEN BLENDERS 35

STERILE WATER FOR MECH. VENTILATION AND INHALATION..... 37

MOUNTS & BRACKETS FOR HUMIDIFIERS 39

WILAFLOW ELITE NIV INFANT CPAP 42

WILAFLOW ELITE NIV INFANT CPAP ACCESSORIES 43

WILABINO NEO POSITIONING AID..... 47

MEDICAL GAS HOSES 48

ACTIVE RESPIRATORY GAS HUMIDIFIERS

AIRcon Gen2 Humidifier



When a patient is requiring mechanical ventilation over a longer period of time, it is imperative to take measures to compensate the losses of heat and moisture in order to avoid the aforesaid complications. Mechanical ventilation with cold and dry respiratory gases causes within short time mucus on the respiratory epithelium to become more viscous, impairing the functionality of the cilia.

AIRcon Gen2 provide a complete respiratory gas humidifying set in combination with its specially designed humidifier chambers and heated wire breathing sets.

It is suitable for mechanically ventilated patients requiring individual therapy with high performance and safety. The humidifier AIRcon Gen2 combines modern technology and innovative design in a high performance unit. Mechanically ventilated patients are provided with optimally conditioned respiratory gas.

Available in selected markets only! Please contact WILAméd for more details.

Consists of:

- Heater Base / Controller
- Reusable Temperature Probe
- Reusable Heater Wire Adapter

Item	Part number	QTY	GTIN
AIRcon Gen2 humidifier 220 - 240 V Version	101200	1	4055146013182
AIRcon Gen2 humidifier 110 - 120 V Version	101201	1	4055146013199

AIRniva Controller



In case of noninvasive ventilation (e.g. respiratory masks) a continuous positive flow is administered (e.g. CPAP). Thereby caused increased oral breathing may elicit undesirable accompanying symptoms. In the long run the upper respiratory tracts dry out through the permanent supply with cool respiratory gases of positive pressure.

The consequences are painfully inflamed nasal and oral mucous membranes as well as blockages of the air passages and congestion of secretion in the respiratory apparatus. Particularly leakages at the respiratory mask may promote drying out of the nasal mucous membranes.

AIRniva provides a complete noninvasive respiratory humidification set in combination with its specially designed humidifier chambers and breathing tube sets.

It is suitable for mechanically, noninvasive ventilated patients requiring individual therapy with high performance and safety. The humidifier AIRniva combines modern technology and innovative design in a high performance unit. Mechanically, noninvasive ventilated patients are provided with optimally conditioned respiratory gas.

Available in selected markets only! Please contact WILAméd for more details.

- Consists of:
- Heater Base / Controller
 - Reusable Heater Wire Adapter INSP

Item	Part number	QTY	GTIN
AIRniva humidifier 220 - 240 V Version	100500	1	4055146012314
AIRniva humidifier 110 - 120 V Version	100505	1	4055146013175

Dual Heater Wire Adaptor INSP + EXP.



Dual heating wire adapter for WILAmEd dual limb sets.
Only for use on the AIRcon humidifier!

Item	Part number	QTY	GTIN
Dual Heater Wire Adaptor INSP + EXP.	100929	1	4055146117705

Single Heater Wire Adaptor INSP only



Heating wire adapter only for inspiratory heated WILAmEd limb sets.
Suitable for use on AIRcon and AIRniva.

Item	Part number	QTY	GTIN
Single Heater Wire Adaptor INSP only	100942	1	4055146800157

Temperature Probe




Suitable for use on AIRcon and AIRniva.

Item	Part number	QTY	GTIN
Temperature Probe 160 cm	100910	1	4055146110560
Temperature Probe 180 cm	100909	1	4055146117712

DISPOSABLE DOUBLE LIMB BREATHING SETS - HEATED


Double-limb set, NEO



Double limb set for AIRcon, heated (i+e), A-adapter, with and without autofill chamber (120cm, 10mm Ø)
Connecting tube 60cm

Item	Part number	QTY	GTIN
BTS1154A, including Auto-Fill humidification chamber	270745	10	4055146032411
BTS1181A, without Auto-Fill humidification chamber	270754	10	4055146030066


Double-limb set, NEO



Double limb set, heated (i+e), A-adapter, with and without autofill chamber (120cm, 10mm Ø)
Connecting tube 60cm

Item	Part number	QTY	GTIN
BTS1122A, including Auto-Fill humidification chamber	270705	10	4055146032237
BTS1125A, without Auto-Fill humidification chamber	270706	10	4055146300015


Double-limb set, Adult



Breathing set heated (i+e) with coiled heated wire for AIRcon, double hose, A-adapter, autofill chamber (150cm, 22mm Ø),

Item	Part number	QTY	GTIN
BTS3481A, including Auto-Fill humidification chamber	271760	10	4055146300787

Double-limb set, Pediatric




Breathing set heated (i+e) with coiled heated wire for AIRcon, double hose, A-adapter, autofill chamber (150cm, 15mm Ø),

Item	Part number	QTY	GTIN
BTS2260A, including Auto-Fill humidification chamber	271761	10	4055146300770

DISPOSABLE NEONATAL HFO BREATHING SETS - HEATED


Double-limb smooth-bore set NEO for Dräger VN500



Breathing set heated (i+e), double hose, for Dräger VN500 HFO, A-adapter, autofill chamber (120cm, 10mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS1245A, including Auto-Fill humidification chamber	271163	10	4055146031292


Double-limb smooth-bore set NEO for SLE 4000/5000



Double limb set for SLE4000/5000, heated (i+e), A-adapter, autofill chamber, pressure measuring line, T-piece with temperature port, flow limiter, extension tube 40cm (120cm, 10mm Ø) connection tube 60cm

Item	Part number	QTY	GTIN
BTS1148A, including Auto-Fill humidification chamber	270739	10	4055146032459


Double-limb smooth-bore set NEO for H&L Leoni Plus



Double limb set for LeoniPlus, heated (i+e), A-adapter, autofill chamber, pressure line, bacteria filter, HFO kit, extension tubing for incubator (120cm, 10mm Ø) connection tube 60cm

Item	Part number	QTY	GTIN
BTS1215A, incl. HFO kit and Auto-Fill humidification chamber	270791	10	4055146030141


Double-limb smooth-bore set NEO for Fabian HFO



Breathing set heated (i+e), double limb, for Acutronic Fabian HFO, A-adapter, autofill chamber AIRcon, additional pressure measurement line (160cm, 10mm Ø), connecting tube 60cm

Item	Part number	QTY	GTIN
BTS1347A, Breathing set heated (i+e), for Acutronic	271710	10	4055146300305

Double-limb smooth-bore set NEO for SLE 6000




Double limb breathing set, heated (i+e), for SLE 6000, A- adapter, autofill chamber AIRcon, additional pressure measuring line (160cm, 10mm Ø), connecting tube 60cm

Item	Part number	QTY	GTIN
BTS1321A SLE 6000,10mm, 160cm, heated	271587	10	4055146300732

DISPOSABLE SINGLE LIMB BREATHING SETS - HEATED


Single-limb set NEO O2 Therapy



Single-limb set for O2 therapy, heated (i), A-adapter, with and without autofill chamber Aircon, Neo adapter, safety valve P40 (120cm, 10mm Ø)
Connecting tube 60cm

Item	Part number	QTY	GTIN
BTS1307A, including Auto-Fill humidification chamber	271506	10	4055146030943
BTS1235A, without Auto-Fill humidification chamber	271507	10	4055146300169


Single-limb set NEO CPAP Therapy



Single-limb set for nCPAP, heated (i), A-adapter, with and without autofill chamber, Neo adapter, nCPAP generator, extension tubing for incubator (120cm, 10mm Ø)
Connecting tube 40cm

Item	Part number	QTY	GTIN
BTS1204A, including Auto-Fill humidification chamber	270777	10	4055146030448
BTS1193A, without Auto-Fill humidification chamber	270771	10	4055146030035


Single-limb set, ADULT



Single-limb set for CPAP, heated (i), A-adapter, with and without autofill chamber (22mm Ø)
Connecting tube 60cm

Item	Part number	QTY	GTIN
BTS3234A, Single-limb set incl. Auto-Fill humidification chamber	271125	10	4055146031063
BTS3237A, Single limb set without Auto-Fill humidification chamber	271126	10	4055146031056
BTS3451A, Single-limb set incl. Universal Auto-Fill humid. chamber	271657	10	4055146300084


Single-limb set, PEDIATRIC



Single limb set for CPAP, heated (i), A-adapter, autofill chamber Aircon (150cm, 15mm Ø)
Connecting tube 60cm

Item	Part number	QTY	GTIN
BTS2171A, Single-limb set incl. Auto-Fill humidification chamber	271519	10	4055146030981


Single-limb set, ADULT (AIRniva only)



Single-limb set for CPAP, O2 therapy, heated (i), A-adapter, autofill chamber universal (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3421A - Single-limb set incl. Auto-Fill humidification chamber	271579	10	4055146033104

Single-limb set ADULT Single use for Philips V60




Single-limb set for BiPAP Vision, V 60, heated (i), A-adapter, autofill chamber (22mm Ø), connection tube 60cm with bacteria filter

Item	Part number	QTY	GTIN
BTS3218A, 150 cm, incl. Auto-Fill humidification chamber	270756	10	4055146030059
BTS3191A, 180 cm, without Auto-Fill humidification chamber	271196	10	4055146030615

DISPOSABLE BREATHING SETS - NON-HEATED

Double-limb smooth-bore set




Standard Ventilation Set (Dry)

2 x 180 cm, 22 mm ID for adults (271003), Smooth lumen
2 x 150 cm, 15 mm ID for children (270329), Smooth lumen

Item	Part number	QTY	GTIN
BTS300, Double-limb smooth-bore set, Adult	271003	10	4055146040935
BTS200, Double-limb smooth-bore set, Pediatric, disposable	270329	10	4055146040294

Coaxial Breathing Circuit



The coaxial (limb-in-limb) ventilation set is a space-saving, easy-to-connect limb set that can be used universally in intensive care as well as in anesthesia ventilation.

In addition to its remarkably low weight and low noise level, the heat transfer and insulation within the set are decisive advantages of this design.

Item	Part number	QTY	GTIN
BTS805, Coaxial Breathing Circuit 200 cm, Adult	270359	10	4055146040591
BTS850, Coaxial Breathing Circuit 300 cm, Adult	270867	10	4055146042250

HUMIDIFICATION CHAMBERS

The humidification chamber is an essential component of our heated ventilation sets, so the two are also perfectly matched.

Thus, we ensure optimal respiratory gas conditioning for mechanically ventilated patients, both in clinical and non-clinical settings.

The C200AF Auto-Fill humidification chambers guarantee a constant water level, on the one hand due to our sophisticated float set inside the chamber and on the other hand due to the highly precise infrared water level monitoring of the AIRcon.

Our humidification chambers can be ordered individually or as part of a complete ventilation set (convenience kit).

With our sterile water „WILaqua“ we can offer a closed set completely from one source.

C200AF AIRcon Disposable Auto-Fill Humidification Chamber



Auto-Fill disposable humidification chamber, for all patient groups.

Only suitable for use on AIRcon humidifiers.

Item	Part number	QTY	GTIN
C200AF, disposable Auto-Fill Humidification Chamber for AIRcon	500300	30	4055146750148
C200AFP, disposable Auto-Fill Humidification Chamber for AIRcon	500480	30	4055146750001


C200AF Universal, disposable Humidification Chamber for AIRniva



Auto-Fill disposable humidification chamber for all patient groups.
For use on PMH1000, PMH5000 and AIRniva and suitable third-party humidifiers only.

Item	Part number	QTY	GTIN
C200AF Universal, disposable Humidification Chamber for AIRniva	500380	30	4055146750049

C200R Universal Humidification Chamber



Manually fillable, reusable humidification chamber for AIRcon and AIRniva humidification chamber, autoclavable


Item	Part number	QTY	GTIN
C200R, reusable Humidification Chamber for AIRcon and AIRniva	500350	1	4055146750032

DISPOSABLE BREATHING SETS HOME CARE - HEATED

For use with Resmed Astral Ventilators




Single-limb set, Adult



Single limb set for Astral 100/150, heated (i), A- adapter, autofill chamber, pressure measuring line, control line, nominal flow 60 l/min ≤ 0.2 kPa (150cm, 22mm Ø) connecting tube 60cm

Item	Part number	QTY	GTIN
BTS3391A, Single-limb set, incl. AF Humidification Chamber	271250	10	4055146032800


Double-Limb Set, Adults



Breathing set heated (i+e) for AIRcon, double hose, A-adapter, autofill chamber (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3167A, Double-limb set incl. AF Humidification Chamber	270678	10	4055146031889


Single-limb set, Pediatric



Single-limb set for Astral 100/150, heated (i), A-adapter, autofill chamber, nominal flow rate 50 l/min ≤ 0.2kPa (150cm, 15mm Ø), connecting tube 60cm

Item	Part number	QTY	GTIN
BTS2200A, Single-limb set, Pediatric, incl. AF Humidification Chamber	271252	10	4055146042267

Double-limb set, Pediatric




Double limb set, heated (i+e), A-adapter, autofill chamber, water trap (e), expiratory side 200cm, small breathing volumes (150cm, 15mm Ø), connecting tube 60cm

Item	Part number	QTY	GTIN
BTS6038A, Double-limb set, Pediatric, incl. AF Humidification Chamber	271508	10	4055146030929

For use with Breas Vivo 50 and Vivo 60 Ventilators




Single-limb set, Adult



Breathing set heated (i), single tube, for Vivo 50/60, A-adapter, autofill chamber (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3238A, Single-limb set, Adult, incl. AF Humidification Chamber	271128	10	4055146031049


Double-limb set, Adult



Breathing set heated (i+e) for AIRcon, double hose, A-adapter, autofill chamber (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3481A, Double-limb set, Adult, incl. AF Humidification Chamber	271760	10	4055146300787


Single-limb set, Pediatric



Breathing set heated (i), single tube, for Vivo 50, Monnal T40, A-adapter, Aircon Autofill chamber (150cm, 15mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS2125A, Single-limb set, Pediatric, incl. AF Humidification Chamber	271132	10	4055146031049

Double-limb set, Pediatric




Breathing set heated (i+e) for AIRcon, double hose, A-adapter, autofill chamber (150cm, 15mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS2260A, Double-limb set, Pediatric, incl. AF Humidification Chamber	271761	10	4055146300770

For use with Trilogy 100 Ventilators




Single-limb set, Adult, AV



Single limb set for Trilogy 100, heated (i), A-adapter, autofill chamber Aircon (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3256A, Single-limb set, Adult, AV, incl. AF Humidification Chamber	271143	10	4055146031414


Single-limb set, Adult, PV



Single limb set for Trilogy 100, Vivo 50/60, heated (i), A-adapter, autofill chamber (150cm, 22mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS3268A, Single-limb set, Adult, PV, incl. AF Humidification Chamber	271156	10	4055146031452


Single-limb set, Pediatric, AV



Single limb set for Trilogy 100, heated (i), A-adapter, Aircon Autofill chamber (150cm, 15mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS2140A, Single-limb set, Pediatric, AV, incl. AF Humidification Chamber	271154	10	4055146031339

Single-limb set, Pediatric, PV




Leakage set, heated (i), A-adapter, autofill chamber (150cm, 15mm Ø), connecting tube 60cm

Item	Part number	QTY	GTIN
BTS2142A, Single-limb set, Pediatric, PV, incl. AF Humidification Chamber	271157	10	4055146031315

For use with Trilogy EVO Ventilators




Double-limb set, Pediatric



Double limb set TEV, heated (i+e), A-adapter, pressure measuring line, autofill chamber (150cm, 15mm Ø), connection tube 60cm

Item	Part number	QTY	GTIN
BTS6078A, Double-limb set, Pediatric, incl. AF Humidification Chamber	271746	10	4055146499528

Double-limb set, Adult



Double limb set TEV, heated (i+e), A-adapter, pressure measuring line, autofill chamber (150cm, 22mm Ø), connection tube 60cm


Item	Part number	QTY	GTIN
BTS8121A, Double-limb set, Adult, incl. AF Humidification Chamber	271747	10	4055146300695

DISPOSABLE BREATHING SETS HOME CARE - NON-HEATED

For use with Resmed Astral Ventilators




Single-limb set, Adult



Single-limb set for Astral 100/150, pressure measuring line, control line, nominal flow 55 l/min ≤ 0.2 kPa (180cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS3390, Single-limb set, Adult	271249	10	4055146042168


Double-Limb Set, Adults



Double limb set (180cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS300 - Double-Limb Set, Adults	271003	10	4055146040935


Single-limb set, Pediatric



Single-limb set for Astral 100/150, pressure measuring line, control line, nominal flow 55 l/min ≤ 0.2 kPa (150cm, 15mm Ø)

Item	Part number	QTY	GTIN
BTS2199 - Single-limb set, Pediatric	271251	10	4055146042267

Double-limb set, Pediatric




Double-limb set (150cm, 15mm Ø)

Item	Part number	QTY	GTIN
BTS200 - Double-limb set, Pediatric	270329	10	4055146040294

For use with Breas Vivo 50 and Vivo 60 Ventilators




Single-limb set, Adult



Breathing set unheated, single tube, for Breas Vivo 50 (180cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS3172, Single-limb set, Adult	270685	10	4055146041413


Double-Limb Set, Adults



Double-limb set (180cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS300 - Double-Limb Set, Adults	271003	10	4055146040935


Single-limb set, Pediatric



Breathing set unheated, single tube, for Breas Vivo 50, Vivo 60 (150cm, 15mm Ø)

Item	Part number	QTY	GTIN
BTS2127, Single-limb set, Pediatric	271136	10	4055146031018

Double-limb set, Pediatric



Double-limb set (150cm, 15mm Ø)

Item	Part number	QTY	GTIN
BTS200 - Double-limb set, Pediatric	270329	10	4055146040294

For use with Cough Assisting Devices



Single-limb set, incl. Bacteria Filter



Single limb set for Comfort Cough, Cough Assist, bacteria filter (90cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS810, Single-limb set, incl. Bacteria Filter	270347	10	4055146040645

Single-limb system, incl. Mask



Single limb set for ComfortCough, 22 MM ID,with bacteria filter and mask (90cm, 22mm Ø)

Item	Part number	QTY	GTIN
BTS253, Single-limb system, incl. Mask für Children	270348	10	4055146040638
BTS349, Single-limb set, incl. Mask for Adults	270349	10	4055146040621
BTS296, Single-limb set XL, incl. Bacteria Filter and Mask for Adults	271192	10	4055146040249

DISPOSABLE BREATHING TUBES HOME CARE - NON-HEATED

The use of safe materials, advanced manufacturing technology and demand- and patient-oriented production have made WILAméd an internationally recognized manufacturer and supplier of breathing tubes. According to our corporate philosophy, customer satisfaction and reliability have top priority. We demonstrate this to our customers on a daily basis, both at home and abroad.

Requirements, suggestions and criticism are always welcome, as this information helps us to adapt our products to changing wishes and needs.

We do not use solutions or chemicals! This makes our breathing hoses not only odorless, but also of high quality and absolutely environmentally friendly.

CPAP Tube



Item	Part number	QTY	GTIN
BT3101, CPAP Tube WILAbasic, gray, 180 cm (6ft.)	2701661	10	4055146020050
BT2100, CPAP Tube WILAlite, gray, 180 cm (6ft.)	2704461	10	4055146020401
Tube Sleeve NORDKAP, gray	270158	10	4055146086384

DISPOSABLE PATIENT INTERFACES

Our extensive range of patient interfaces for non-invasive ventilation therapy is full of ideas and innovations and offers optimal comfort for the patient, which clearly contributes to the success of the therapy.

Ventura Total Full Face Masks



WILAmEd’s Ventura Total Face Masks are designed and manufactured with patient comfort in mind. Our VENTURA breathing masks are less intimidating than other full face masks because they are open-cut, allowing for better treatment outcomes for anxious or claustrophobic patients.

Item	Part number	QTY	GTIN
Ventura TFM NV/SBE (Blue), Size S	201112	1	4055146600085
Ventura TFM NV/SBE (Blue), Size M	201111	1	4055146600092
Ventura TFM NV/SBE (Blue), Size L	201110	1	4055146600016
Ventura TFM CPAP/AAV (Clear), Size S	201113	1	4055146600023
Ventura TFM CPAP/AAV (Clear), Size M	201114	1	4055146600030
Ventura TFM CPAP/AAV (Clear), Size L	201115	1	4055146600047
Ventura TFM NV/AAV (Yellow), Size S	201107	1	4055146600139
Ventura TFM NV/AAV (Yellow), Size M	201108	1	4055146600146
Ventura TFM NV/AAV (Yellow), Size L	201109	1	4055146600153

WILAhiflow Mini Plus Nasal Cannulas



WILAmEd offers a wide range of nasal cannulae for use in nasal high-flow oxygen therapy: from low flow rates when ventilating the smallest patients to very high flow rates for adults.

Item	Part number	QTY	GTIN
WILAhiflow Mini Plus Nasal Cannula (NC11S) Size S	30076715	15	948538362518
WILAhiflow Mini Plus Nasal Cannula (NC11M) Size M	30076815	15	948538362501
WILAhiflow Mini Plus Nasal Cannula (NC11L) Size L	30076915	15	948538362495
WILAhiflow Mini Plus NoSlip-Pads	30077100	100	-

High-Flow Nasal Cannulas



Specifically designed for Nasal High Flow Therapy, the nasal cannula allows delivery of flow rates up to 80 liters per minute. Its compact, flexible design, together with its light weight and low noise, provides maximum patient comfort.

Item	Part number	QTY	GTIN
High-Flow Nasal Cannula, Size S	20123310	10	16970531333735
High-Flow Nasal Cannula, Size M	20123410	10	16970531333728
High-Flow Nasal Cannula, Size L	20123510	10	16970531333711

LMFlo2 High-Flow Nasal Cannulas



WILAméd offers a wide range of nasal cannulae for use in nasal high-flow oxygen therapy: from low flow rates when ventilating the smallest patients to very high flow rates for adults.

Item	Part number	QTY	GTIN
LMFlo2 High Flow Nasal Cannula, Size S (Adult)	201230	10	4719876361202
LMFlo2 High Flow Nasal Cannula, Size M (Adult)	201231	10	4719876361196
LMFlo2 High Flow Nasal Cannula, Size L (Adult)	201232	10	4719876360403

Catheter Mount, Double Swivel



Tube extension with angled double rotary connector and suction bronchoscopy opening. Fits all common endotracheal tubes, laryngeal masks and laryngeal tubes.

Item	Part number	QTY	GTIN
Catheter Mount, Double Swivel, 20 cm, 22M/15F-22F	2701031	31	4055146450109
Catheter Mount, Double Swivel, 20 cm, 22M/15F-15M	2705001	31	4055146450079
Catheter Mount, Double Swivel, 15 cm, 22M/15F-22F	2705011	31	4055146450086
Catheter Mount, Double Swivel, 15 cm, 22M/15F-15M	2701021	31	4055146450093

HEAD GEAR FOR PATIENT INTERFACES - HOME-CARE

Cap: DALLAS



Item	Part number	QTY	GTIN
Cap DALLAS, Navy Blue, Size M	300012	10	4055146070086
Cap DALLAS, Navy Blue, Size L	300013	10	4055146070154

Cap: NEW YORK




Item	Part number	QTY	GTIN
Cap NEW YORK, White, Size S	300005	10	4055146071052
Cap NEW YORK, White, Size M	300010	10	4055146071090

Head Strap: RIO, MILANO and BARCELONA



Item	Part number	QTY	GTIN
Head Strap RIO, Navy Blue, One-Size	300450	10	4055146070192
Head Strap MILANO, Turquoise, 20mm	201001	10	4055146070581
Head Strap BARCELONA, Turquoise, 20mm	203001	10	4055146070369

Chin Strap: TORONTO and OTTAWA




Item	Part number	QTY	GTIN
Chin Strap TORONTO, White, One-Size	300315	10	4055146071045
Chin Strap OTTAWA, Navy Blue, One-Size	300316	10	4055146071038

MECHANICAL AIR / OXYGEN BLENDERS

WILAmEd/Biomed Devices oxygen mixers for neonatal and pediatric intensive care units with flowmeters from 3 to 15 LPM, for low and high flow oxygen therapy.

The WILAmEd/Biomed Devices High-Flow Oxygen Blenders are specially designed for Nasal High-Flow Therapy in combination with active respiratory gas humidification (WILAmEd AIRcon). The blender has a flowmeter up to 70 liters per minute and is therefore the perfect solution in Nasal High-Flow Therapy for adults.

Biomed Blender



Item	Part number	QTY	GTIN
Biomed NEOblend, including Flowmeter 15 LPM	100356	1	00813841020478
Biomed High-Flow Blender incl. Flowmeter 70 LPM	100357	1	00813841020461

Biomed Blender with 2 Flowmeters



Item	Part number	QTY	GTIN
Biomed Blender incl. 70 LPM and 15 LPM Flowmeter	100361	1	00813841020959

DEHAS Quality Mix



Item	Part number	QTY	GTIN
DEHAS Quality Mix AIR/O2 Blender-Set 85 LPM Flow meter and Silencer	100374	1	4251411700683

STERILE WATER FOR MECH. VENTILATION AND INHALATION

WILAqua Sterile Water



WILAqua sterile water is produced in state-of-the-art manufacturing facilities and has been specially developed to meet the demanding hygienic requirements of invasively and non-invasively ventilated patients. It is distilled sterile water.

Item	Part number	QTY	GTIN
WILAqua Sterile Water for Inhalation 1.000 ml PE Bottle	500186	10	05200702400596
WILAqua Sterile Water for Inhalation 1.000 ml PE Bottle	5001862	660	05200702400596

ALLEGRA Sterile Water-Pack



ALLEGRA is a closed, disposable sterile water set for simple, hygienic humidification during inhalation and long-term oxygen therapy.

Item	Part number	QTY	GTIN
ALLEGRA Sterile Water-Pack for Inhalation, 450 ml PE Bottle	5001866	25	9555768100413
ALLEGRA Sterile Water-Pack for Inhalation, 450 ml PE Bottle	5001867	1100	9555768100413

MOUNTS & BRACKETS FOR HUMIDIFIERS

IV Pole / C-Clamp



Our standard c-clamp is designed to mount horizontally or vertically on round surface (e.g. I.V. poles), with a maximum diameter of 40 mm.

Item	Part number	QTY	GTIN
IV Pole / C-Clamp	550301	1	4055146087244

Mounting Bracket for Norm Rail



Our standard mounting unit is designed, for example, for mounting the AIRcon humidifier on standard rails.

Item	Part number	QTY	GTIN
Mounting Bracket for Norm Rail	550226	1	4055146087848

Mounting Bracket Kit for Norm Rail



Our angled attachment unit is intended, for example, for attaching the AIRcon respiratory gas humidifier to standard rails, including the INTENSA ventilator.

Item	Part number	QTY	GTIN
Mounting Bracket Kit for Norm Rail	550247	1	4055146088562

Telescopic Water Bag Pole



Water bottle holder with telescopic rod and holder for respiratory gas humidifier, for attachment to standard rails in clinical and non-clinical environments.

Item	Part number	QTY	GTIN
Telescopic Water Bag Pole	550220	1	4055146085967

Biomed C-Clamp



Biomed C-clamp is designed for fixing e.g. the oxygen mixer to infusion stands, with a maximum diameter of 40 mm.

Item	Part number	QTY	GTIN
Biomed Blender Pole Mount C-Clamp	100358	1	4055146200742

Heavy Duty Stand IV-Pole



Heavy duty infusion stand made of stainless steel for use as a mobile equipment carrier such as oxygen mixers, respiratory gas humidifiers, etc., including 4 hooks for attaching sterile water, infusion solutions, etc.

Item	Part number	QTY	GTIN
Heavy Duty Stand with five Casters and IV-Pole	550273	1	4055146087749

WILAFLOW ELITE NIV INFANT CPAP

WILAflow Elite NIV Ventilator



The WILAflow Elite is a ventilator that is used for non-invasive ventilation, such as nCPAP and nasal HiFlow O2 therapy, in premature and newborn infants in the hospital.

The use of nCPAP (a continuous positive airway pressure) is now considered the gold standard for respiratory support in preterm infants. The most common indications include respiratory distress syndrome and apneas. Application via binasal prongs or nasal mask is now considered the internationally validated standard.

WILAflow Elite ventilation modes: SNIPPV (Synchronized), NIPPV, nCPAP, HFNC

All accessories and consumables for the WILAflow Elite are available and are described on the following pages.

- Scope of delivery:
- WILAflow Elite ventilator
 - tubing support arm
 - Supply tree for sterile water
 - Signal Box/IEC

Item	Part number	QTY	GTIN
WILAflow Elite NIV Ventilator, incl. Vent Cart	101300	1	4055146013267

WILAFLOW ELITE NIV INFANT CPAP ACCESSORIES

WILAflow Elite Signal Box



Item	Part number	QTY	GTIN
WILAflow Elite Signal Box for Abdominal Sensor	101302	1	4055146800119

WILAflow Elite Body Sensor




Item	Part number	QTY	GTIN
WILAflow Elite Body Sensor VX010 100 CM	101301	5	4055146800171

WILAflow Elite Circuit Support Arm



Item	Part number	QTY	GTIN
WILAflow Elite Circuit Support Arm	101304	1	4055146800133

WILAflow nCPAP Generator



The Fluidic-Flip principle, which has been known and proven for decades when using our CPAP generator, was developed especially for infants and has established itself as the globally recognized gold standard in the non-invasive ventilation of infants, or premature infants.

Our WILAflow nCPAP generator comes with three nasal prongs; 1 x size S, 1 x size M and 1 x size L.

Item	Part number	QTY	GTIN
WILAflow nCPAP Generator, incl. 3 Prongs (S, M, L)	300660	10	4055146083079

WILAflow Nasal Prongs



Item	Part number	QTY	GTIN
WILAflow Nasal Prong, Size XS	300732	10	5060195764700
WILAflow Nasal Prong, Size S	300712	10	4055146640111
WILAflow Nasal Prong, Size M	300713	10	4055146640159
WILAflow Nasal Prong, Size L	300714	10	4055146640104

WILAflow CPAP Nasal Mask



Item	Part number	QTY	GTIN
WILAflow CPAP Nasal Mask, Size S	300710	10	4055146200261
WILAflow CPAP Nasal Mask, Size M	300711	10	4055146200438
WILAflow CPAP Nasal Mask, Size L	300708	10	4055146200513
WILAflow CPAP Nasal Mask, Size XL	300709	10	4055146200285

WILAbonnets



For use in combination with our WILAflow nCPAP generator, our particularly soft, flexible and skin-friendly WILAbonnets are used for optimal fixation of the nCPAP generator during nCPAP therapy.

Item	Part number	QTY	GTIN
WILAbonnet, White, Size 000, 18 - 20 cm	300720	10	4055146200216
WILAbonnet, Grey, Size 00, 20 - 22 cm	300721	10	4055146200223
WILAbonnet, Pink, Size 0, 22 - 24 cm	300722	10	4055146200230
WILAbonnet, Brown, Size 1, 24 - 26 cm	300723	10	4055146200247
WILAbonnet, Yellow, Size 2, 26 - 28 cm	300724	10	4055146200254
WILAbonnet, Blue, Size 3, 28 - 30 cm	300725	10	4055146200612
WILAbonnet, Lite Orange, Size 4, 30 - 32 cm	300726	10	4055146200674
WILAbonnet, Green, Size 5, 32 - 34 cm	300727	10	4055146200537
WILAbonnet, Red, Size 6, 34 - 36 cm	300728	10	4055146200148
WILAbonnet, Orange, Size 7, 36 - 38 cm	300729	10	4055146200049
WILAbonnet, Turquoise, Size 8, 38 - 40 cm	300730	10	4055146200032
WILAbonnet, Dark Blue, Size 9, 40 - 42 cm	300731	10	4055146200155

WILABINO NEO POSITIONING AID

WILAbino NEO



WILAbino NEO - safe, developmental and durable head positioning set as a stabilization aid for premature and newborn babies. WILAbino NEO consists of various components that can be used in different variations. The basic product consists of two components: A support ring and an outer protective cover made of fabric.

Item	Part number	QTY	GTIN
WILAbino NEO Positioning TG PUR GEL, Size S	201302	6	4260764840571
WILAbino NEO Positioning TG PUR GEL, Size M	201301	6	4260764840564
WILAbino NEO Positioning TG PUR GEL, Size L	201300	6	4260764840557
WILAbino NEO HSR Single Patient Use Cover, Size S	201305	50	4260764840663
WILAbino NEO HSR Single Patient Use Cover, Size M	201304	50	4260764840656
WILAbino NEO HSR Single Patient Use Cover, Size L	201303	50	4260764840649

MEDICAL GAS HOSES

O2 and AIR



ZV hoses for WILAmEd ventilators for the compressed air and oxygen in ISO32 black and white.

Our ZV hoses convince with durable and uncompromising reliability and offer absolute safety for man and machine.

Item	Part number	QTY	GTIN
O2 Medical Gas Hose, NIST/DIN, white, 3 m	101577	1	4055146800027
AIR Medical Gas Hose, NIST/DIN, black, 3 m	100666	1	4055146083703

Space for your notes

[illegible]